

## Tilla Care Product Catalogue

### General Sales information

UriCap Female Devices and Extension Tubes SX-L can be ordered separately in bulk or packaged in kits and boxed according to the product list below.



US patent No.:  
11,666,474

FDA Registered Class I  
Medical Device


### Components:

The UriCap Female is assembled from two parts, the UriCap Female and the Extension Tube.

Products in the catalogue and in the kits has been tested with the UriCap Female for optimum performance and functionality.



### UriCap Female Device specifications

### Extension Tube SX specifications



Name	UriCap Female Device	Name	Extension Tube SX-L
<b>PN</b>	FPR02020L	<b>PN</b>	FPR00004
<b>Material</b>	Medical grade silicone	<b>Material</b>	Tube- Medical grade silicone Connector- Medical grade polypropylene
<b>Size packed</b>	6*4*1 in. 15*10*2.5 cm.	<b>Size packed</b>	5*5.5*0.6 in. 25*7.5*1.5 cm.
<b>Weight packed</b>	0.24 oz 7g	<b>Tube Length</b>	11 in. 28 cm
<b>Shelf life</b>	3 years	<b>Weight packed</b>	0.33 oz 9.6 gr
<b>Storage conditions</b>	50-86°F 10-30°C	<b>Shelf life</b>	3 years
<b>US Reimbursement code (HCPCS)</b>	A4328	<b>Storage conditions</b>	50-86°F 10-30°C
<b>Comments</b>	Single Use device. certified for maximum 24 hrs. continuous use	<b>US Reimbursement code (HCPCS)</b>	A4331
		<b>Comments</b>	The Extension Tube should be replaced with the urine bag

### Product List:

Image	Part #	Description	Quantities	Weight	Size
	FPR02020L	UriCap Female Device	Bulk is minimum order of 1000 pcs	247 oz 7 kg.	23*20.5*16 in. 58*52*41 cm.
	FPR00004	Extension Tube SX-L	Bulk is minimum order of 250 pcs	77 oz 2.5 kg.	19*14.5*11 in. 49*37*28 cm.

#### Local distributor:

#### Manufactured by:

**TillaCare Ltd.**



2 Etgar st., 390321 Tirat Carmel,  
Israel

E-mail: [sales@tillacare.com](mailto:sales@tillacare.com)







#### For more information:

[www.tillacare.com/training](http://www.tillacare.com/training)

Image	Part #	Description	Quantities	Weight	Size
	HRD00006	Securement Device	50	11.1 oz 315 gr.	7.3*7.3*3.7 in. 18.5*18.5*9.5 cm
	HRD00004	Lubricating Jelly	144	21.3 oz 604 gr.	3.3*5.5*5.1 in. 8.5*14*13 cm.
	HRD00007	Alcohol Prep Pads	200	6.2 oz 175 gr.	4.3*4.1*2.1 in.
	KIT00010	UriCap Female Demonstration Kit	-	11.8 oz 335 gr.	8.8*3.8*2.2 in 24*14*9 cm
	HRD00008	Urine Bag Tilla 2L	35	70.55 oz 2 kg	13*9*8 in 33*23*20.5 cm.

#### Kit Product List:

Image	Part #	Description	Kit Items	Quantities	Size
	KIT00011	Starter Kit	<p>1pc <b>FPR02020L</b> UriCap Female device</p> <p>1pc <b>FPR00004</b> Extension Tube SX-L</p> <p>1pc <b>HRD00004</b> Sterile Lubricating Jelly</p> <p>1pc <b>HRD00006</b> Securement device</p> <p>1pc <b>HRD00007</b> Alcohol Prep Pad</p>	<p>35 kits 35.27 oz 1 kg</p>	<p>10.5*8*7 in. 27*21*19 cm</p>
	KIT00013-30	Base Kit	<p>1pc <b>FPR02020L</b> UriCap Female device</p> <p>1pc <b>FPR00004</b> Extension Tube SX-L</p> <p>1pc <b>HRD00004</b> Sterile Lubricating Jelly</p>	<p>30 kits 24.7 oz 700 g</p>	<p>10.5*8*7 in. 27*21*19 cm</p>

	KIT00012	Change Kit	<p>1pc <b>FPR02020L</b> UriCap Female device</p> <p>1pc <b>HRD00004</b> Sterile Lubricating Jelly - US</p>	<p>70 kits 28.2 oz 800g</p>	<p>10.5*8*7 in. 27*21*19 cm</p>
	KIT00020-30	Base Kit+	<p>1pc <b>FPR02020L</b> UriCap Female device Assembled to:</p> <p>1pc <b>FPR00004</b> Extension Tube SX-L</p> <p>1pc <b>HRD00004</b> Sterile Lubricating Jelly</p>	<p>30 kits 24.7 oz 700 g</p>	<p>10.5*8*7 in. 27*21*19 cm</p>